



ELSEVIER

# Rapidly inducible, genetically targeted inactivation of neural and synaptic activity *in vivo*

Dougal Tervo<sup>1</sup> and Alla Y Karpova<sup>2</sup>

Inducible and reversible perturbation of the activity of selected neurons *in vivo* is critical to understanding the dynamics of brain circuits. Several genetically encoded systems for rapid inducible neuronal silencing have been developed in the past few years offering an arsenal of tools for *in vivo* experiments. Some systems are based on ion-channels or pumps, others on G protein coupled receptors, and yet others on modified presynaptic proteins. Inducers range from light to small molecules to peptides. This diversity results in differences in the various parameters that may determine the applicability of each tool to a particular biological question. Although further development would be beneficial, the current silencing tool kit already provides the ability to make specific perturbations of circuit function in behaving animals.

## Addresses

<sup>1</sup> Watson School of Biological Sciences, Cold Spring Harbor Laboratory, 1 Bungtown Rd, Cold Spring Harbor, NY 11724, United States

<sup>2</sup> Howard Hughes Medical Institute, Janelia Farm Research Campus, 19700 Helix Drive, Ashburn, VA 20147, United States

Corresponding author: Karpova, Alla Y ([karpovaa@janelia.hhmi.org](mailto:karpovaa@janelia.hhmi.org))

**Current Opinion in Neurobiology** 2007, 17:581–586

This review comes from a themed issue on  
New technologies  
Edited by Mark Schnitzer and Atsushi Miyawaki

Available online 28th November 2007

0959-4388/\$ – see front matter

© 2007 Elsevier Ltd. All rights reserved.

DOI 10.1016/j.conb.2007.10.002

## Introduction

Underlying the function and dysfunction of the brain are complex neuronal circuits, consisting of diverse classes of specific neuronal populations. A major goal in systems neuroscience is to relate the activity in specific neuronal circuits to physiological processes, behavioral responses, and disease states. This is a challenging task especially in the case of the highly complex mammalian brain and requires tools that allow rapidly inducible and reversible modulation of specific parts of neural circuits *in vitro* and *in vivo*; a need that has been widely recognized for some time [1–4]. The past decade has seen an upsurge in the development and optimization of genetically encodable tools for perturbing neuronal activity and several systems have been described [5<sup>•</sup>,6,7,8<sup>••</sup>,9<sup>•</sup>,10,11<sup>••</sup>]. More importantly, initial *in vivo* applications of these tools have been described in the past two years opening the exciting new field of highly specific perturbation of circuit

function in behaving animals [8<sup>••</sup>,12<sup>••</sup>,13<sup>••</sup>,19<sup>••</sup>]. As these systems differ in the approach to neuronal silencing, mechanism of triggering, and kinetics of induction and reversal, each may be best applied in different experimental contexts. In this review, we discuss the most promising rapidly inducible silencing systems with the focus on how the different properties of each system make it particularly suitable for certain types of biological questions while highlighting the need for further characterization and optimization of their various features.

## Approaches to rapid silencing


### Allatostatin receptor-allatostatin system

Callaway and colleagues pioneered the use of the *Drosophila* allatostatin receptor (AlstR) [14] for inducible inactivation of neuronal activity (Figure 1A) [6]. This system is based upon the expression of a G protein coupled receptor that is selectively activated by the insect peptide hormone allatostatin (AL) but is insensitive to endogenous mammalian peptides [14]. Following application of allatostatin, AlstR indirectly causes opening of G protein inward rectifier K<sup>+</sup> (GIRK) channels leading to a hyperpolarizing potassium current, which in turn suppresses action potential generation. In dissociated cultures, a decrease in excitability is observed within 10–20 min of allatostatin application. Similarly, in slices of the spinal cord from mice expressing AlstR in a subpopulation of spinal cord inhibitory neurons, reversible suppression of excitability was observed in the presence of allatostatin (10 nM) with a time course of approximately 20 min for both induction and reversal [15<sup>•</sup>]. Tan *et al.* used this system successfully in rats, ferrets, and monkeys to selectively suppress neuronal activity in genetically defined visual cortical and thalamic neurons [12<sup>••</sup>]. Local application of allatostatin (100 nM) to the relevant neuronal tissue (within ~1 mm) was necessary because of restricted diffusion of the peptide, and its effect on neuronal excitability was inducible and reversible with the similar time course of 10–20 min.

### The ivermectin-sensitive chloride channel

A second system that reduces neuronal activity upon ligand application is the ivermectin (IVM)-sensitive chloride channel consisting of *Caenorhabditis elegans* GluCl  $\alpha$  and  $\beta$  subunits (Figure 1B) [7,13<sup>••</sup>,16–18]. Expression of GluCl $\alpha/\beta$  using Sindbis virus in dissociated hippocampal neurons allowed inducible and potent inactivation of neuronal activity within seconds of ivermectin (5 nM) application. The unbinding of ivermectin from the channel appeared to be very slow as reversal *in vitro* required 8 h. Although the wild type GluCl $\alpha/\beta$  channel was also

Figure 1



Schematics of five rapidly inducible and reversible, genetically encoded neuronal inactivation systems described here. **(A)** Allatostatin receptor (AlstR; green) couples the binding of allatostatin (red spheres) to the activation of a G protein (white). This then leads to the opening of a G protein coupled inward rectifier K<sup>+</sup> channel (blue). **(B)** The two subunit GluCl $\alpha$ / $\beta$  receptor (blue, purple) is opened by ivermectin (black circle). The resulting Cl<sup>-</sup> current suppresses neuronal activity. **(C)** Chloride current through the five subunit GABA<sub>A</sub> channel is enhanced by binding of zolpidem (black circle) to the interface between  $\gamma$ 2 and  $\alpha$ 1 subunits leading to suppression of excitability. **(D)** Two Molecular Systems for the Inactivation of Synaptic Transmission (MISTs). The one-component system is based on crosslinking of modified VAMP (light blue). The two-component system is based on mislocalization of modified Synaptophysin (Sph; purple) non-specifically along the plasma membrane using the transmembrane domain of Syntaxin (StxTM). **(E)** Halorhodopsin (NpHR) pumps chloride into the cell when opened by yellow light.

sensitive to glutamate, this side-effect was subsequently reduced [18]. Lerchner *et al.* demonstrated that this system can work *in vivo* by using it to perturb the function of striatal circuits in rats [13<sup>••</sup>]. By expressing GluCl $\alpha$ / $\beta$  unilaterally in the striatum through co-infection with two adeno-associated viruses encoding the individual subunits, they demonstrated that intraperitoneal administration of 5–10 mg of ivermectin per kilogram of body weight led to a perturbation of amphetamine-induced rotational behavior. The observed effect was maximal at 12–24 h post injection and reversed by 4 days post injection. It is important to emphasize that unlike in the *in vivo* study with the AlstR/AL system, the readout of the inactivation in this case was not a direct measure of the electrical activity of the targeted neuronal population, but rather behavioral output, making it difficult to compare the time course of silencing between the two systems. It is not unlikely that the change in behavioral output occurs on a slower time scale than the change in the activity of targeted neurons.

#### The zolpidem-sensitive GABA<sub>A</sub> channel

The most recently developed system takes advantage of the specificity of the GABA<sub>A</sub> channel binding to a positive allosteric modulator, zolpidem [19<sup>••</sup>]. Single amino acid substitution in one of the subunits of the GABA<sub>A</sub> channel ( $\gamma$ 2) leads to abolition of zolpidem sensitivity. Wulff *et al.* used this result to create a strategy for cell type specific enhancement of GABA-ergic currents by first creating zolpidem-insensitive mice using a mutant  $\gamma$ 2 subunit and then reintroducing the wild type subunit to genetically

defined neuronal populations to restore the sensitivity to the drug to select neuronal populations [19<sup>••</sup>]. In acute cerebellar slices from mice expressing the zolpidem-sensitive version of the GABA<sub>A</sub> receptor selectively in Purkinje cells, enhancement of miniature inhibitory postsynaptic currents (mIPSC) in Purkinje cells was observed within minutes of zolpidem (1  $\mu$ M) delivery. Furthermore, intraperitoneal administration of the inducer in these transgenic mice led to impairment of rotarod performance within 5 min. Neither the extent of block of the excitatory input onto Purkinje cells by zolpidem administration nor the time course of reversal was addressed.

#### MISTs

A completely different approach to silencing is based on blocking neurotransmission by perturbing the function of synaptic proteins; an approach originally exemplified by *shibire*, a temperature-sensitive mutant of *Drosophila* dynamin [20]. A different embodiment of this strategy that has proven effective in the mammalian system is Molecular Systems for the Inactivation of Synaptic Transmission (MISTs) [8<sup>••</sup>]. MISTs interfere with the synaptic transmission by inducing the mislocalization or immobilization of modified presynaptic proteins using chemically inducible dimerization (CID) with small cell-permeable chemicals, which are analogs of rapamycin [21]. Two systems have been introduced: VAMP/Syb MIST is based on oligomerization of synaptic vesicle protein VAMP/Synaptobrevin, whereas Sph-StxTM MIST is based on inducible crosslinking of synaptic vesicle protein Synaptophysin to the plasma membrane through

the transmembrane domain of Syntaxin (Figure 1C). It was demonstrated that in the presence of the relevant dimerizers (100 nM homodimerizer and 250 nM heterodimerizer) the two MISTs induced a reversible block of synaptic transmission in cultured neurons and brain slices within 20–30 min. It was further shown that the VAMP/Syb MIST system can work *in vivo*. In transgenic mice expressing VAMP/Syb MIST selectively in Purkinje neurons, intraventricular administration of 0.5 nmol of the homodimerizer impaired learning of and performance in the rotarod behavioural task. A decrease in behavioral performance was observed on the first session following dimerizer administration (~12–18 h post administration); the performance recovered 2 days later. The time course of the actual silencing of the targeted neurons *in vivo* was not determined.

### Rhodopsin/halorhodopsin

All four of the abovementioned approaches depend upon the application of an exogenous chemical ligand. An alternative approach is controlling neuronal activity by manipulating neuronal excitability using light-gated membrane channels or pumps (Figure 1D). One example of this is rhodopsin 4 (RO4) that partially suppresses action potential generation by hyperpolarizing neurons in a light-dependent manner within milliseconds [9<sup>\*</sup>]. Another potent system is based on the light-activated chloride pump halorhodopsin (NpHR), which, in combination with the cation permeable light-sensitive channelrhodopsin 2 (ChR2), permits bi-directional control of action potential firing with millisecond resolution [10,11<sup>••</sup>]. The applicability of halorhodopsin to *in vivo* manipulation of circuit activity has been verified by demonstration of the control of movement in *C. elegans* [11<sup>••</sup>]. The rhodopsin and NpHR systems hold a great promise for rapid *in vivo* manipulation of neural activity.

### What to consider when choosing the silencing approach

The relevant biological question and experimental paradigm will define several criteria that will guide the choice of the silencing method. A non-exhaustive list of the important issues to consider includes the relevant time frame for the induction and reversal of silencing; the nature of the neuronal population to be targeted for

inactivation; the desired expression strategy and the constraints on inducer administration.


### Time course of inactivation and its reversal

Different experimental paradigms will impose distinct constraints on the time course of induction, maintenance, and reversal of silencing (Figure 2). For instance, some experiments in sensory physiology rely on *in vivo* whole-cell recording of subthreshold membrane activity [22,23]. The limited recording time available in such experiments requires that the circuit perturbation be induced and reversed within minutes. Thus, of the existing methods halorhodopsin might be best suited in this context if silencing needs to be reversible, whereas GABA<sub>A</sub> or GluCl $\alpha/\beta$  can also be used when reversal is not important. Experiments examining the functional structure and dynamics of cortical circuits might employ imaging of neuronal activity using functional indicators [24,25]. Although a number of genetically encodable indicators have been described [26,27], their inability to follow single action potentials has limited their applicability up to now [28]. Thus, indicator dyes that are bulk-loaded into tissue are still employed by most studies [25,29,30]. The time course of such experiments is limited by the retention half-life for these synthetic dyes and usually ends up being on the orders of several hours. In such cases, circuit perturbation needs to be achieved within tens of minutes for longitudinal experiments; thus, all of the abovementioned systems should be applicable. In contrast to network dynamics questions, studies exploring the contribution of particular circuit components to animal behavior will often tolerate a longer induction time course for activity perturbation. All of the systems described in this review should be useful in such studies. In cases when long-term silencing is required, however, the exogenous conductances in channel-based approaches may perturb the energy budget of neurons possibly making MISTs a better choice.

### The nature of circuit components to be silenced

It is also important to consider which circuit components are to be targeted. One issue is whether the neuronal population is widely distributed or localized to a small region of the brain. As the allatostatin receptor and halorhodopsin require localized delivery of the inducer

Figure 2



Time course of the induction and reversal of inactivation for the four neuronal silencing systems.

(allatostatin or light), they may be best suited for the latter case. By contrast, the GABA<sub>A</sub>, GluCl $\alpha/\beta$ , and the MIST systems should be applicable in both. In some cases the limited spread of the induction by AlstR and light will be a great advantage if the required degree of cell type specificity cannot be achieved by genetic means. Another issue is how large the cells are physically and/or electronically. When a high level of expression for exogenous channels is difficult to achieve genetically, channel density might decrease with the size of the targeted cell. Thus, small cells, for example, interneurons, could be easily silenced by channel-based methods whereas large cells, for example, Purkinje cells, might prove difficult to hyperpolarize and may be better targeted for inactivation using MISTs. On the contrary, if neurons in the circuit are electrically coupled, channel-based methods would be advantageous as MISTs selectively perturb communication through chemical synapses. Suppression of excitability through the zolpidem/GABA<sub>A</sub> system depends on the interaction of the drug with the  $\alpha$ - $\gamma$  subunit interface of the channel. Zolpidem is most effective at the  $\alpha$ 1- $\gamma$ 2 subunit containing complexes with a significantly lower affinity for the  $\alpha$ 2 and  $\alpha$ 3 containing channels and no affinity for the  $\alpha$ 4- $\alpha$ 6 subunits [31,32]. As the zolpidem/GABA<sub>A</sub> systems rely on the endogenous alpha subunit expression, silencing of neurons in regions with little or no  $\alpha$ 1 subunit expression (like amygdala, ventral striatum, and others [33]) may prove challenging. Finally, at present, MISTs represent the only approach for perturbing individual projections within a circuit. Since MIST-dependent inactivation of neuronal activity relies on the perturbation of synapses, localized delivery of the dimerizer should allow targeting of individual projections. By contrast, channel-based interference with action potential generation will lead to inactivation of the entire set of projections from the targeted neuronal population. It is important to mention, however, that if good axonal targeting of channel-based systems can be achieved, localized inhibition of action potential propagation leading to silencing of specific projections may become possible for these systems as well.

### Expression strategy

Two strategies for introducing silencing systems into animals have been pursued: viral expression with stereotactic delivery [12<sup>••</sup>,13<sup>••</sup>] and transgenesis [8<sup>••</sup>]. For both strategies the choice between one and two-component systems is important. Whereas two-component systems (GluCl $\alpha/\beta$  and Sph-Stx<sup>TM</sup> MIST) will give a greater combinatorial control over the population of neurons to be silenced, the accompanying requirement for co-expression of both components can be technically difficult in viral and genetic systems. One-component systems are easier to express and allow for much simpler breeding schemes when used in transgenic mice, which is an important issue when many animals are needed for experiments.

### Delivery of the inducer

Another factor to consider is the possible restriction on the way in which the inducer is delivered. For AlstR and halorhodopsin, the delivery has to be rather localized. Allatostatin needs to be applied within  $\sim$ 1 mm of the targeted neurons. If cells in the cortex are being silenced, the peptide can be applied to the surface of the brain; however, it needs to be injected with a glass pipette if deeper brain areas are targeted [12<sup>••</sup>]. The light source for the activation of halorhodopsin may need to be close to the targeted neurons as well given the power required for stimulation and the light-scattering and absorption properties of the brain tissue and is often done with an optical fiber for stimulating deep in the brain. Although both glass pipettes for AL delivery and optical fibers can be quite small in diameter ( $\sim$ 30  $\mu$ m for the pipette [12<sup>••</sup>] and  $\sim$ 100  $\mu$ m for the fiber [34]) and thus may not cause much damage, care may need to be taken in some cases to ensure that an important fiber tract has not been disrupted. Dimerizers that are used for the induction of MIST-dependent silencing do not impose such strict restrictions on delivery. Karpova *et al.* administered the homodimerizer used with VAMP/Syb MIST into the lateral ventricle, far removed from the terminals of the Purkinje cells [8<sup>••</sup>]. Further studies will be needed to determine if this homodimerizer penetrates the blood-brain barrier (BBB) sufficiently well that it can be administered peripherally. Ivermectin and zolpidem are so far the only inducers discussed here that have been shown to work in the brain following peripheral delivery [13<sup>••</sup>]. The heterodimerizer used with the Sph-Stx<sup>TM</sup> MIST is known to cross the blood-brain barrier (Dr Victor Rivera, personal communication), and thus should be able to work when administered systemically as well.

For many experiments, the stress to animals following surgical delivery of the inducer may prove to be a non-issue. In others, a single surgery implanting the optical fiber or a miniature osmotic pump days in advance of the actual behavioral testing will circumvent the problem of stress. Yet in some cases, the behavioral studies may be exquisitely sensitive to any stress the animal experiences. Under those circumstances, the less invasive the route of the inducer delivery is, the more advantageous is the silencing system. In the ideal case, the chemical inducer would be delivered with food or water and the inducer would not only need to cross the blood-brain barrier but also be sufficiently water soluble. This can be accomplished for the tetracycline-inducible gene expression, which has been useful for the slower induction of neuronal silencing using tetanus toxin light chain [35] and could be used to make the elegant system based on tethered toxins [5<sup>•</sup>] inducible as well. However, none of the inducers discussed here are sufficiently water soluble for this route of delivery. Short of that, intraperitoneal injection would be the best drug delivery method when stress is a factor.

## Current challenges

Although great progress has been made developing genetically encoded systems for rapidly inducible and reversible silencing, none of the systems described here is ideal. The limitations of the existing systems can be partly overcome by a careful choice of the silencing system (or systems) given the biological question and experimental context. However, further characterization of (1) the pharmacokinetics of the inducers, (2) the time course of the actual neuronal silencing rather than of behavioral perturbation, and (3) the dynamic range of the effect at different expression levels of each system would help in making a proper choice. In addition, it is already clear that some optimization of the existing systems would be beneficial. For instance, it is known that ivermectin can cause toxicity at high concentrations [36]. Although no toxicity was observed in the study by Lerchner *et al.*, the concentration of ivermectin necessary to achieve behavioral perturbation was only twofold to fourfold smaller than what has been reported to cause side-effects [6]. Furthermore, ivermectin affects other endogenous murine channels including GABA channels [37] and P2X4 purinergic receptors [38]. It will be important to address these potentially negative side-effects in the future. Similarly, the heterodimerizer used with the Sph-StxTM MIST system binds to the endogenous FKBP protein and is thus likely to have side-effects at some concentrations. Finding ways to modify these chemicals to alleviate these non-specific effects while retaining their functionality would be a major advance. Improving their blood–brain barrier permeability may be helpful as well.

## Conclusions

With the recent developments in the field of molecular tools for perturbing neural activity, highly specific manipulations of circuit function *in vivo* are becoming possible. Although some challenges remain, we are finally well under way to go beyond classical lesion studies in understanding the neural basis of behavior.

## Acknowledgements

The authors thank Winfried Denk, Thomas Stüdhof, and Karel Svoboda for helpful comments on the manuscript. DGRT was supported, partly, by the Howard Hughes Medical Institute predoctoral fellowship and George A and Majorie H Anderson Fellowship. AYK was supported, partly, by the Burroughs Wellcome Foundation and Howard Hughes Medical Institute (HHMI). Development of the MIST silencing system was supported by HHMI and NIH (R01MH070052 to Karel Svoboda).

## References and recommended reading

Papers of particular interest, published within the annual period of review, have been highlighted as:

- of special interest
  - of outstanding interest
1. Shah EM, Jay DG: **Methods for ablating neurons.** *Curr Opin Neurobiol* 1993, **3**:738-742.
  2. Marek KW, Davis GW: **Controlling the active properties of excitable cells.** *Curr Opin Neurobiol* 2003, **13**:607-611.
  3. Miesenbock G: **Genetic methods for illuminating the function of neural circuits.** *Curr Opin Neurobiol* 2004, **14**:395-402.
  4. Callaway EM: **A molecular and genetic arsenal for systems neuroscience.** *Trends Neurosci* 2005, **28**:196-201.
  5. Ibanez-Tallon I, Wen H, Miwa JM, Xing J, Tekinay AB, Ono F, Brehm P, Heintz N: **Tethering naturally occurring peptide toxins for cell-autonomous modulation of ion channels and receptors *in vivo*.** *Neuron* 2004, **43**:305-311.
  6. Lechner HA, Lein ES, Callaway EM: **A genetic method for selective and quickly reversible silencing of mammalian neurons.** *J Neurosci* 2002, **22**:5287-5290.
  7. Slimko EM, McKinney S, Anderson DJ, Davidson N, Lester HA: **Selective electrical silencing of mammalian neurons *in vitro* by the use of invertebrate ligand-gated chloride channels.** *J Neurosci* 2002, **22**:7373-7379.
  8. Karpova AY, Tervo DG, Gray NW, Svoboda K: **Rapid and reversible chemical inactivation of synaptic transmission in genetically targeted neurons.** *Neuron* 2005, **48**:727-735. Provides *in vitro* and *in vivo* characterization of MISTs.
  9. Li X, Gutierrez DV, Hanson MG, Han J, Mark MD, Chiel H, Hegemann P, Landmesser LT, Herlitze S: **Fast noninvasive activation and inhibition of neural and network activity by vertebrate rhodopsin and green algae channelrhodopsin.** *Proc Natl Acad Sci U S A* 2005, **102**:17816-17821. Characterizes rhodopsin 4-dependent silencing *in vitro* and in chick embryos.
  10. Han X, Boyden ES: **Multiple-Color Optical Activation, Silencing, and Desynchronization of Neural Activity, With Single-Spike Temporal Resolution.** 2007: e299.
  11. Zhang F, Wang L-P, Brauner M, Liewald JF, Kay K, Watzke N, Wood PG, Bamberg E, Nagel G, Gottschalk A *et al.*: **Multimodal fast optical interrogation of neural circuitry.** *Nature* 2007, **446**:633-639. Introduces the use of the chloride pump halorhodopsin for rapid suppression of neuronal activity and demonstrates its use in *C elegans*.
  12. Tan EM, Yamaguchi Y, Horwitz GD, Gosgnach S, Lein ES, Goulding M, Albright TD, Callaway EM: **Selective and quickly reversible inactivation of mammalian neurons *in vivo* using the *Drosophila* allatostatin receptor.** *Neuron* 2006, **51**:157-170. Extends the use of allatostatin receptor to *in vivo* applications in the thalamus and cortex of rats, ferrets, and monkeys.
  13. Lerchner W, Xiao C, Nashmi R, Slimko EM, van Trigt L, Lester HA, Anderson DJ: **Reversible silencing of neuronal excitability in behaving mice by a genetically targeted, ivermectin-gated *Cl(-)* channel.** *Neuron* 2007, **54**:35-49. *In vivo* application of the GluCl system to silencing of the striatal circuits.
  14. Birgül N, Weise C, Kreienkamp HJ, Richter D: **Reverse physiology in drosophila: identification of a novel allatostatin-like neuropeptide and its cognate receptor structurally related to the mammalian somatostatin/galanin/opioid receptor family.** *EMBO J* 1999, **18**:5892-5900.
  15. Gosgnach S, Lanuza GM, Butt SJB, Saueressig H, Zhang Y, Velasquez T, Riethmacher D, Callaway EM, Kiehn O, Goulding M: **V1 spinal neurons regulate the speed of vertebrate locomotor outputs.** *Nature* 2006, **440**:215-219.
  16. Cully DF, Vassilatis DK, Liu KK, Pareiss PS, Van der Ploeg LH, Schaeffer JM, Arena JP: **Cloning of an ivermectin-sensitive glutamate-gated chloride channel from *Caenorhabditis elegans*.** *Nature* 1994, **371**:707-711.
  17. Li P, Slimko EM, Lester HA: **Selective elimination of glutamate activation and introduction of fluorescent proteins into a *Caenorhabditis elegans* chloride channel.** *FEBS Lett* 2002, **528**:77-82.
  18. Slimko EM, Lester HA: **Codon optimization of *Caenorhabditis elegans* GluCl ion channel genes for mammalian cells dramatically improves expression levels.** *J Neurosci Methods* 2003, **124**:75-81.
  19. Wulff P, Goetz T, Leppa E, Linden AM, Renzi M, Swinny JD, Vekovischeva OY, Sieghart W, Somogyi P, Korpi ER *et al.*: **From synapse to behavior: rapid modulation of defined neuronal**

- types with engineered GABA(A) receptors. *Nat Neurosci* 2007, **10**:923-929.  
Describes the use of zolpidem-sensitive GABA<sub>A</sub> channel to enhance IPSCs onto targeted cells.**
20. Kitamoto T: **Targeted expression of temperature-sensitive dynamin to study neural mechanisms of complex behavior in Drosophila.** *J Neurogenet* 2002, **16**:205-228.
  21. Spencer DM, Wandless TJ, Schreiber SL, Crabtree GR: **Controlling signal transduction with synthetic ligands.** *Science* 1993, **262**:1019-1024.
  22. Margrie TW, Brecht M, Sakmann B: **In vivo, low-resistance, whole-cell recordings from neurons in the anaesthetized and awake mammalian brain.** *Pflugers Arch* 2002, **444**:491-498.
  23. Lee AK, Manns ID, Sakmann B, Brecht M: **Whole-cell recordings in freely moving rats.** *Neuron* 2006, **51**:399-407.
  24. Ohki K, Chung S, Ch'ng YH, Kara P, Reid RC: **Functional imaging with cellular resolution reveals precise micro-architecture in visual cortex.** *Nature* 2005, **433**:597-603.
  25. Ohki K, Chung S, Kara P, Hubener M, Bonhoeffer T, Reid RC: **Highly ordered arrangement of single neurons in orientation pinwheels.** *Nature* 2006, **442**:925-928.
  26. Mank M, Reiff DF, Heim N, Friedrich MW, Borst A, Griesbeck O: **A FRET-based calcium biosensor with fast signal kinetics and high fluorescence change.** *Biophys J* 2005, **90**:1790-1796.
  27. Miyawaki A, Griesbeck O, Heim R, Tsien RY: **Dynamic and quantitative Ca<sup>2+</sup> measurements using improved cameleons.** *Proc Natl Acad Sci U S A* 1999, **96**:2135-2140.
  28. Pologruto TA, Yasuda R, Svoboda K: **Monitoring neural activity and [Ca<sup>2+</sup>] with genetically encoded Ca<sup>2+</sup> indicators.** *J Neurosci* 2004, **24**:9572-9579.
  29. Diez-Garcia J, Akemann W, Knopfel T: **In vivo calcium imaging from genetically specified target cells in mouse cerebellum.** *Neuroimage* 2007, **34**:859-869.
  30. Sullivan MR, Nimmerjahn A, Sarkisov DV, Helmchen F, Wang SS: **In vivo calcium imaging of circuit activity in cerebellar cortex.** *J Neurophysiol* 2005, **94**:1636-1644.
  31. Pritchett DB, Seeburg PH: **Gamma-aminobutyric acidA receptor alpha 5-subunit creates novel type II benzodiazepine receptor pharmacology.** *J Neurochem* 1990, **54**:1802-1804.
  32. Wafford KA, Thompson SA, Thomas D, Sikela J, Wilcox AS, Whiting PJ: **Functional characterization of human gamma-aminobutyric acidA receptors containing the alpha 4 subunit.** *Mol Pharmacol* 1996, **50**:670-678.
  33. Laurie DJ, Wisden W, Seeburg PH: **The distribution of thirteen GABAA receptor subunit mRNAs in the rat brain. III. Embryonic and postnatal development.** *J Neurosci* 1992, **12**:4151-4172.
  34. Mehta AD, Jung JC, Flusberg BA, Schnitzer MJ: **Fiber optic in vivo imaging in the mammalian nervous system.** *Curr Opin Neurobiol* 2004, **14**:617-628.
  35. Yamamoto M, Wada N, Kitabatake Y, Watanabe D, Anzai M, Yokoyama M, Teranishi Y, Nakanishi S: **Reversible suppression of glutamatergic neurotransmission of cerebellar granule cells in vivo by genetically manipulated expression of tetanus neurotoxin light chain.** *J Neurosci* 2003, **23**:6759-6767.
  36. Dadarkar SS, Deore MD, Gatne MM: **Comparative evaluation of acute toxicity of ivermectin by two methods after single subcutaneous administration in rats.** *Regul Toxicol Pharmacol* 2007, **47**:257-260.
  37. Adelsberger H, Lepier A, Dudel J: **Activation of rat recombinant alpha(1)beta(2)gamma(2S) GABA(A) receptor by the insecticide ivermectin.** *Eur J Pharmacol* 2000, **394**:163-170.
  38. Sim JA, Chaumont S, Jo J, Ulmann L, Young MT, Cho K, Buell G, North RA, Rassendren F: **Altered hippocampal synaptic potentiation in P2X4 knock-out mice.** *J Neurosci* 2006, **26**:9006-9009.